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The screenshot shows a control quality interface. At the top, there are three buttons: 'PageUp', 'PageDown', and 'Open'. The 'PageUp' and 'PageDown' buttons are highlighted with a red dashed border. Below the buttons, there is a text field containing the lot number 'LOT: 99512-28.2-0141-0223-14-880809-0'. Below the text field is a table with three columns: 'Item Name', 'Result', and 'Unit'. The table contains the following data:

Item Name	Result	Unit
TP	69.9	g/L
ALB	40.0	g/L
TBIL	10.0	umol/L
ALT	30	U/L
AST	20	U/L
GGT	20	U/L
UREA	5.00	mmol/L

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The screenshot shows a mobile application interface for quality control. At the top, there are two buttons: 'Back' on the left and 'Search' on the right. Below these buttons is a list of two items, each with a label and a date, and an 'edit' button with a pencil icon. The 'fromDate' is '01/01/2000' and the 'toDate' is '01/01/2021'. Red dashed boxes highlight the 'Search' button and the 'edit' buttons for both items.

fromDate:	01/01/2000	edit
toDate:	01/01/2021	edit

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The screenshot shows a user interface for quality control. At the top, there are two buttons: "No" on the left and "Yes" on the right. Below these buttons is a large text area containing the question: "Determine reagent disc Liver and Renal Function?". This question is enclosed in a red dashed rectangular box. The background of the interface is light gray.

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The screenshot displays the MMDMP Medical Data Management Platform interface. The 'Setting' menu is highlighted with a red dashed box and a circled '1'. A secondary window titled 'Setting' is open, showing various configuration options. The 'Quality Control Parameters' option is highlighted with a red dashed box and a circled '2'. The main interface includes sections for 'Inspection Info', 'Patient info', and 'Result'.

Inspection Info

Analysis Date: 2021/08/23
 Patient Area: Out-patient
 ID:
 Name:
 Sex: Male
 Age:
 Department: ICU
 Bed No.:
 Sample Type: Serum
 Submitter: 送检人1
 Submit Time: 2021/08/23 14:09
 Operator: 检验人1
 Reviewer: 审核人1
 Diagnosis: 传染病
 Remark: 建议出院

Patient info

Name	Sex	Age	ID	Time	Device number

Result

Item	Result	Ranges	Unit

Setting

Basic Setting
 Software Info | Language | Baud Rate Setting | **Quality Control Parameters**

Report Setting
 Hospital | Add Items | **Quality Control Parameters**

Test Info Setting
 Patient Area | Sex | Department | Sample Type
 Submitter | Operator | Reviewer | Diagnosis
 Remark | Submit Time

Legend:
 ** Indicates that the result is beyond the dynamic range
 Yellow: Add test item from other device
 Blue: Below lower limit | Red: Above upper limit | Green: Out of value

Buttons: Select all, Display All, Preview(V), Print(F2), Combine, LIS upload, Delete, Add Item

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Quality Control Parameters

Disc Type: Glucose and Lipid Panel Lyophilized Kit [5]

QC Batch Number: p0616001

Search

Save

Name	Sex	Age	ID	Time	Device number	Item	Result	Ranges	Unit
GLU						0.00	0.00		mmol/L
TG						0.00	0.00		mmol/L
CHOL						0.00	0.00		mmol/L
HDL-C						0.00	0.00		mmol/L
S-GSP						0.00	0.00		mmol/L

** Indicates that the result is beyond the dynamic range

■ Add test item from other device
■ Below lower limit ■ Above upper limit ■ Out of value

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The screenshot shows the MNCHIP Medical Data Management Platform interface. The 'Result' menu item in the top navigation bar is highlighted with a red dashed box. A secondary window titled 'Result' is open, showing three buttons: 'Patient', 'Quality Control Result', and 'Quality Control Chart'. The 'Quality Control Chart' button is also highlighted with a red dashed box. The main interface includes sections for 'Inspection Info', 'Patient info', and 'Result'.

Inspection Info

Analysis Date: 2021/08/23
 Patient Area: Out-patient
 ID:
 Name:
 Sex: Male
 Age:
 Department: ICU
 Bed No.
 Sample Type: Serum
 Submitter: 送检人1
 Submit Time: 2021/08/23 14:09
 Operator: 检验人1
 Reviewer: 审核人1
 Diagnosis: 没病
 Remark: 建议出院

Patient info

Name	Sex	Age	ID	Time	Device number

Result

Item	Result	Ranges	Unit

Buttons: Patient, Quality Control Result, Quality Control Chart

Legend: ** indicates that the result is beyond the dynamic range
 Add test item from other device
 Below lower limit (purple), Above upper limit (red), Out of value (green)

Bottom bar: Select all, Display All, Preview(V), Print(F2), Combine, LIS upload, Delete

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The screenshot displays the MNCHIP Medical Data Management Platform interface. The main window is titled "Quality control Result" and contains several sections:

- Inspection Info:** A form on the left with fields for Analysis Date (2021/08/23), Patient Area (Out-patient), ID, Name, Sex (Male), Age, Department (ICU), Bed No., Sample Type (Serum), Submitter (送检人1), Submit Time (2021/08/23 14:09), Operator (检验人1), Reviewer (审核人1), Diagnosis (检测), and Remark (建议出检). A "Save" button is at the bottom.
- Patient info:** A table with columns: Name, Sex, Age, ID, Time, Device number.
- Result:** A table with columns: Item, Result, Ranges, Unit.
- Search Filters (highlighted in red dashed box):**
 - Starting time: 2021/08/23
 - End Time: 2021/08/23
 - Device number: All
 - QC Batch Number: All
 - Search button
- Legend:**
 - ** Indicates that the result is beyond the dynamic range
 - Yellow: Add test item from other device
 - Purple: Below lower limit
 - Red: Above upper limit
 - Green: Out of value
 - Add Item button
- Bottom Bar:** Select all, Display All, Preview(V), Print(F2), Combine, LIS upload, Delete.

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The screenshot shows the MNChip Medical Data Management Platform (v1.0.4.27.4) interface. The main window is titled "Quality Control Chart" and contains the following elements:

- Inspection Info Panel:**
 - Analysis Date: 2021/09/07
 - Patient Area: Out-patient
 - ID: [Empty]
 - Name: [Empty]
 - Sex: Male
 - Age: [Empty]
 - Department: Anorectal department
 - Bed No.: [Empty]
 - Sample Type: Serum
 - Submitter: 送检人1
 - Submit Time: 2021/09/07 16:18
 - Operator: 检验人1
 - Reviewer: 审核人1
 - Diagnosis: 结肠癌
 - Remark: 建议出院
- QC Data Entry Form:**
 - QC Date: 2019/09/07 (To: 2021/09/07)
 - Device Number: 41460
 - Disc Type: Glucose and Lipid Panel Lyophilized Kit [5]
 - QC Batch Number: p0616001
 - Item: GLU
 - Search button
- QC Chart:**

Batch	Batch Date	Result
1	06-16 6.00	+3SD
2	06-16 5.00	+3SD
3	06-16 8.00	+3SD
4	06-16 15.0	+3SD

The chart shows a target line at approximately 0 and control limits at +3SD and -3SD. The data points for batches 1, 2, and 3 are near the target, while batch 4 is significantly higher.
- Bottom Panel:**
 - Buttons: Select all, Display All, Preview(V), Print(F2), Combine, Delete
 - Legend:
 - Yellow: Add test item from other device
 - Blue: Below lower limit
 - Red: Above upper limit
 - Green: Out of value
 - Add Item button

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MNCHIP Technologies Co., Ltd